

# Setting up NobelProcera® Full Contour Zirconia Implant Crown in DTX Studio™ Lab



Final product: NobelProcera Full Contour Zirconia Implant Crown with ASC – base level\*

\*Base is placed at the time of surgery


## DTX Studio™ Lab



### Setting up a case for a Nobel Biocare N1™ Implant

**Note** The following steps show how to order a NobelProcera® Full Contour Zirconia Implant Crown in DTX Studio™ Lab; for details on ordering other restoration types, refer to the appendix

**Note** Please make sure to update to DTX Studio™ Lab 1.12.3 or later; update to the latest version at [DTX Studio™ Go bit.ly/dtxstudiogo](http://bit.ly/dtxstudiogo)



1. Select the **+ sign** in the top left corner to start a new case

The screenshot shows the DTX Studio Lab interface. A drop-down menu is open, with 'Design restoration' highlighted in green. The menu also includes 'Produce model from scan', 'Scan and send model', and 'Import case'. In the background, there is a table of cases with columns for case name, patient name, and status. The cases listed are:

Case Name	Patient Name	Status
pro namio	Dr Nobel	Created
carro	Herminio Fernandes	Scanned
don chon	Herminio Fernandes	Created


2. Select *Design restoration* in the drop-down menu

The screenshot shows the 'Add case' form in DTX Studio Lab. The form is divided into two main sections: 'Case information' and 'Clinician information'. Both sections are highlighted with an orange box. The 'Case information' section includes fields for Patient (a dropdown menu), Case ID (a text input), and Due date (a date picker set to 29 June 2023). The 'Clinician information' section includes fields for Clinician name (a dropdown menu), ZIP code (a text input), and Dental lab or clinic (a dropdown menu). At the bottom of the form, there is a 'Next' button with a right-pointing arrow, which is also highlighted with an orange box.

3. Fill out *Case information* and *Clinician information*

4. Click *Next*


# Selecting details for the case



5. Select **Tooth number** to be restored (example: tooth 14)


6. Under Screw retained Restorations, select **Full contour implant crown**

(For cutback design, define cutback in Design stage)



7. Under Connection type, select **N1 Base Tri**

(This indicates that it is a Nobel Biocare N1™ Implant – trioval conical connection – with the Nobel Biocare N1™ Base placed at the time of surgery)




8. Diameter is automatically set to **RP** (since this selection is in the posterior)

**Note** Diameter is automatically set to **RP** for teeth in the posterior (1-3, 14-19 and 30-32); all other teeth have an option to select **NP**

9. Select **Shade** of your choice

10. Choose **Design flow**  
 - Start from library  
 - Start from diagnostics


## Reviewing and completing the case



11. Review system rendering (example: Nobel Biocare N1™ Base is on the Nobel Biocare N1™ trioval conical connection)

The image of the base indicates the restoration will be designed on top of the Nobel Biocare N1™ Base

The rectangular outline indicates it can utilize an angulated screw channel



notes...




12. Click **Finish** to return to the DTX Studio™ Lab cockpit

You can upload STL files in the DTX Studio™ Lab cockpit to start designing your implant crown

# Appendix/reference guide

## Nobel Biocare N1™ crown restoration selections in DTX Studio™ Lab



**1** **Base** section is for Universal Abutment restorations (TiBase style restorations)

**Note** Titanium abutment blanks are not yet available in DTX Studio™ Lab

