



Surgical Handpiece Zygoma

Product overview

nobelbiocare.com/drillunits



The solution for zygomatic implant application, offering ease of access to treatment sites, safe and effortless handling, and local cooling

Designed for easy access to treatment site

The angulation of the handpiece is designed for optimized access in complex procedures in the zygomatic arch area.



Safe handling

The 20:1 handpiece with contra-angle chuck system combines safe handling with high concentricity.



Effortless handling

The ergonomic design of the handpiece ensures fatigue-free handling.



Local cooling

The external cooling spray emerges at the tip of the instrument and cools the surgical site.





Technical data

Transmission ratio	20:1
Motor connection	ISO 3964
Rotary instruments	Ø 2.35 mm
Maximum drive speed	50'000 rpm
Maximum torque at rotary instrument	70 Ncm
Maximum torque at rotary instrument using rotary instrument with hexagon	105 Ncm

Surgical Handpiece Zygoma

NB30095000	Surgical Handpiece Zygoma SZ-75, 20:1
------------	---------------------------------------

ORDER ONLINE

Order the complete range of our implants and prefabricated prosthetics 24 hours a day through the Nobel Biocare online store.

nobelbiocare.com/store

ORDER BY PHONE

Call our customer service team, or contact your sales representative.

LIFETIME WARRANTY

On all Nobel Biocare implants including prefabricated prosthetic components. For further information visit:

nobelbiocare.com/warranty

nobelbiocare.com/drillunits



87484 GB 2104 Printed in the EU © Nobel Biocare Services AG, 2020. All rights reserved. Distributed by: Nobel Biocare. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Please contact the local Nobel Biocare sales office for current product assortment and availability. Disclaimer: For prescription use only. Caution: Federal (United States) law or the law in your jurisdiction may restrict this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions. The Zygoma handpiece is a device that is manufactured by W&H and is distributed by Nobel Biocare Services AG, Switzerland.